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For professional dental use.

DESCRIPTION

N'Durance® Cristal is a **light-cured, low-shrinkage, high conversion, radio-opaque, esthetic, nanohybrid dental composite material. Its improved translucent properties make it useful for both anterior and posterior restorations.**

N'Durance® Cristal uses a state-of-the-art nano-dimer conversion technology® and nanofiller system.

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This technology provides low volumetric shrinkage for low stress on the bond interface. It exhibits high polymerization conversion with short curing times for improved physical and mechanical properties and less leaching of any residual monomer for better biocompatibility.

N'Durance® Cristal offers superior esthetics and mechanical properties for simplified restorative procedures.

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N'Durance® Cristal is available in the following VITA®* shades: A1, A2, A3, A3.5, B1, B2, C2. It also features an Incisal shade. N'Durance® Cristal provides the dentist with a material that enables the restoration to blend and match the natural appearance of the tooth.

COMPOSITION

EBPDMA (Ethoxylated Bis-Phenol Dimethacrylate)
UDMA (Urethane Dimethacrylate)
DDCDMA (Dimer Dicarbamate Dimethacrylate)

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The state-of-the-art nanofiller system contains 65% filler by volume (10nm – 500nm): ytterbium fluoride mono-dispersed nano-sized particle combined with regular silica and barium glass.

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INDICATIONS

N'Durance® Cristal is indicated for direct esthetic restoration primarily in the anterior region.

This low shrinkage nanohybrid composite can be used in all direct dental restorations, Class I, II, III, IV, and V, and diastema closure. Its body shades match the VITA®* Shade Guide and can be placed using a single layer technique. A final layer of N'Durance® Cristal Incisal shade will add a translucent vitality to the restoration.

CONTRA-INDICATIONS

N'Durance® Cristal is contra-indicated for use with patients who have a history of severe allergic reaction to dimethacrylate resins.

DIRECTIONS FOR USE

1. Cavity cleaning

Teeth should be cleaned by scaling and prophylaxis using flour of pumice without fluoride.

2. Shade selection

Shade selection should be performed prior to isolation and/or preparation of the teeth, ideally in broad daylight on a hydrated tooth. Use of a VITA®* Shade Guide is recommended.

3. Tooth isolation

Isolate the teeth. The use of a rubber dam is highly recommended.

4. Cavity preparation

Use conservative cavity preparations for all Class restorations.

Leave no residual material or base from any previous restoration.

Bevel enamel margin for the anterior restorations.

5. Placement of the matrix

Place the appropriate matrix or crown form to assure proper proximal contour of the tooth. Place wedges to produce good gingival adaptation and tooth separation. Burnish the matrix to achieve optimum proximal contact. The matrix may be placed following the adhesive application, if preferred.

6. Pulp protection

For deep cavities, the pulp should be protected by a cavity liner. Biodentine™- bioactive dentin substitute - is recommended both for pulp exposure and to protect the pulp in deep cavities. Calcium hydroxide cement may also be used for covering small pulp exposures. Glass-ionomer lining cements may also be used in deep cavities. **Eugenol-containing cements must be avoided.**

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Light cure composite with Natural Translucency

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Light Cure Composite with Natural Translucency

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7. Tooth conditioning and bonding

N'Durance® Cristal can be used with any bonding agent. Condition the surface if needed. Follow the manufacturer’s instructions for the adhesive application. Cure the bonding resin before placing the composite, in accordance with the manufacturer’s instructions.

Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of N'Durance® Cristal.

8. Placement of N'Durance® Cristal

Dispense N'Durance® Cristal composite, following the directions corresponding to the dispensing system chosen to restore the tooth:

- Syringe

Dispense onto a mixing pad the needed amount of composite. To prevent more of the material from dispensing when dispensing is completed, aim the tip of the syringe upwards and turn the plunger counter-clockwise to release the pressure and immediately recap the syringe.

Protect the dispensed material from direct light exposure with an opaque cover to avoid premature polymerization. Avoid placing under the direct intensity of the operatory light during the delivery process.

Use a suitable placement instrument.

- Single Dose Capsules

Insert the capsule into a dispenser (e.g. Septodont Multi-Use Dispenser). Deliver an increment of the composite directly into the preparation. Use steady pressure. Rotate the single dose capsule’s tip to the proper angle of entrance into the cavity. Maintain pressure on the dispenser handle while removing it from the mouth to avoid loosening of the capsule. To remove the used single dose capsule, be sure that the dispenser plunger is completely pulled back. Place the used capsule in the proper receptacle.

Placement of a compatible flowable liner, such as N'Durance® Dimer Flow, prior to placement of N'Durance® Cristal is optional. If desired, follow manufacturer’s directions for use.

N'Durance® Cristal composite should be placed into the preparation in increments no greater than 2.5 mm. It is important to use thinner layers for darker shades. Place the composite into the cavity using an appropriately designed dental hand instrument and adapt the composite to the walls of the preparation. Light cure layer by layer. Add increments until the cavity preparation is filled. N'Durance® Cristal is a high viscosity composite that will resist sticking to instruments and slumping, yet will have sufficient flow to adapt to the walls of the cavity preparation and allow sculpting of an anatomical form.

As a final step, a layer of N'Durance® Cristal Incisal shade will add a translucent vitality to the restoration.

9. Photopolymerization

Light cure each layer for 30 seconds with a suitable light-curing lamp (recommendation: use a Quartz-tungsten-halogen curing light with a minimum radiance of 500mW/cm² in the spectral range of 470-480 nanometers). Refer to curing light manufacturer’s recommendations for compatibility and curing recommendations. Do not allow any contamination with saliva or rinsing between cured layers. The N'Durance® Cristal restoration should be exposed to an additional light curing to all surfaces following removal of the matrix to ensure complete polymerization of the material.

10. Finishing and polishing

Begin finishing immediately after the restoration has been completely cured. Remove gross excess to adjust restoration to proper margins, contours, contacts and occlusion.

Finish using a fine diamond, multifluted carbide finishing bur, or finishing discs and strips.

To achieve a very high luster with N'Durance® Cristal, use rubber cups and points with polishing paste.

After the polishing procedure is completed, the composite surface can be cleaned and the surrounding enamel re-etched with phosphoric acid, washed and dried thoroughly. The bonding agent is applied in a thin layer, dried and cured, ensuring that microscopic surface imperfections are removed.

WARNINGS AND PRECAUTIONS FOR USE

- N'DURANCE CRISTAL contains dimethacrylate resins. Avoid the use of this product on patients with known dimethacrylates allergies. To reduce the risk of allergic response, avoid exposure to uncured resin.
- Do not use eugenol-containing materials for pulp protection since they can retard the curing process.
- Do not allow saliva or water contamination of etched tooth surface or unset

material. Use of a rubber dam is highly recommended to achieve optimal results.

- Dentists and assistants should wear gloves and protective eye wear. Patients also should wear eye protection.
- Any non-sterilizable item that is handled in the delivery of the dental service should be disinfected by standard dental office hygiene procedures.
- This product is intended to be used only as described in these Instructions for Use. Any other use is at the sole responsibility of the practitioner.

STORAGE

- N'Durance® Cristal is designed to be stored at ambient temperature, between 50° - 75°F (10° - 24°C).
- Replace cap on syringe immediately after dispensing.
- Do not store in proximity to eugenol-containing products.
- Do not store the product in intense light, or under wet conditions.
- Do not use after expiration date.

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* VITA® is a registered trademark of VITA Zahnfabrik, Bad Säckingen, Germany.