Septo-etch

Etchant 40%

FOR PROFESSIONAL USE ONLY

Instructions for Use

INDICATIONS FOR USE:

Septo-etch is indicated for initiation of the bonding sequence for bonding restorative materials to enamel and/or dentin.

DESCRIPTION:

Septo-etch is a viscous gel containing 40% phosphoric acid. It has a carefully balanced viscosity to achieve the following goals:

1. To allow easy dispensing from the fine metal syringe tip onto the tooth.
2. To have the gel stay in the location where it is placed.
3. To allow the dispensed gel to quickly adapt to the dentin and enamel surfaces on which it is placed.

The gel is supplied in a syringe for easy delivery. Small bore metal syringe tips placed on the syringe allow dispensing of exact amounts of gel directly onto the surfaces to be etched.

FORMULATION

Septo-etch is a water based 40% phosphoric acid gel.

INSTRUCTIONS:

1. The tooth should be cleaned by scaling and prophylaxis with flour of pumice. Isolate tooth and prepare any necessary cavity form according to normal dental procedures.
2. Apply Septo-etch to the tooth surface where etching of the enamel and dentin is required.
3. Allow etch to remain on the surface for 20 seconds. Agitation of the gel with soft brush bristles may enhance the etching process. Do not rub any hard material against the tooth during the etching process.

4. Rinse the etchant away. Start by using a few drops of water with high volume suction to avoid splashing the gel onto the tissues. Then wash with copious amounts of water using an oil-free air/water spray for several seconds.
5. Do not allow saliva contamination of the etched surface.
6. Follow instructions of the bonding agent or resin to be applied for management of the etched surface after the etching step.
7. Remove the syringe tip after use and recap the syringe with the original cap to prevent drying of the gel.

WARNINGS:

Dentists and assistants should wear gloves and protective eye wear. Patients also should wear eye protection.

Etchant gels are corrosive and can cause tissue damage in contact with skin and eyes. If accidental contact occurs with the skin, flush immediately with copious amounts of water. For eye contact, flush immediately and for a prolonged period with copious amounts of water and consult a physician. Do not use in individuals with allergies to these materials.

Any non-sterilizable item that is handled in the delivery of the dental service, should be disinfected by standard dental office hygiene procedures.

STORAGE:

The system is designed to be stored between 59°-75° F (15°-24° C).

Replace caps immediately after dispensing.

If stored in a refrigerator bring to room temperature prior to use.

Manufactured by SEPTODONT
416 S. Taylor Ave | Louisville, CO 80027
800-872-8305 | www.septodontusa.com
Made in the USA

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