# SeptoBond

## LIGHT CURE ADHESIVE

Total etch bonding system

FOR PROFESSIONAL USE ONLY

INSTRUCTIONS FOR USE

### INDICATIONS:

SeptoBond is indicated for bonding to both etched enamel and etched dentin in the restoration of teeth with direct composite resin, compomers, or resin modified glass ionomers.

### DESCRIPTION:

SeptoBond is a light-curing adhesive resin system designed to bond to enamel and dentin. Used along with Septo-etch etching gel, it can provide a quick and easy way to bond composite restorations to teeth.

### FORMULATION:

SeptoBond is a dimethacrylate based monomer system and photosensitive catalyst in an acetone solvent.

### INSTRUCTIONS:

- The tooth should be cleaned by scaling and prophylaxis with flour of pumice. Isolate tooth with a rubber dam and prepare cavity form according to normal dental procedures.
- Base as necessary using a thin layer of calcium hydroxide material on pulp exposures. The calcium hydroxide can be covered with a thin layer of glass ionomer base cement. Avoid use of any eugenol containing cements.
- Apply a 40% phosphoric acid gel such as Septo-etch to the tooth surface for 20 seconds. Enamel can be etched longer by carefully applying the etching agent first to the enamel margins of the preparation. (Teeth with high fluoride content and primary teeth may require longer enamel etch times.) Then apply to dentin making sure to cover all cut dentinal surfaces.
- IMPORTANT: Rinse the dentin quickly with water within 15 -20 seconds of application of the acid etchant.
- Continue to wash with an oil-free air/water spray for 10 – 15 seconds.
- Allow water to stay on dentin until ready to apply SeptoBond. The dentin should be moist, not dehydrated, and without any pooling of water.
- Dispense 1-2 drops of SeptoBond into the mixing well just prior to application. Replace cap immediately after dispensing. Use a brush to carry material to the tooth.
- Just before applying SeptoBond to the tooth, quickly blow with oil-free, dry air for a second. Then apply the bonding agent thoroughly to all walls of the preparation being careful not to rub hard against etched enamel.
- Allow the solvent to evaporate for 30 seconds with a very gentle stream of oil-free, dry air.
- Confirm that all surfaces to be bonded appear glossy. If not, reapply SeptoBond and dry to achieve glossy appearance.
- Using a dental curing lamp with an output around 470-480 nanometers and a minimum of 400mW/cm<sup>2</sup> intensity, light cure for 20 seconds. The bonding agent is compatible with both the conventional quartz-tungsten-halogen curing light and the LED curing light.
- Place the composite in layers of about 2.5 mm or less. Refer to manufacturer's instructions for its placement, curing and finishing.

### WARNINGS:

- Dentists and assistants should wear gloves and protective eye wear. Patients also should wear eye protection.
- Etchant gels are caustic and can cause tissue damage in contact with skin or eyes. If accidental contact with skin occurs, flush immediately with copious amounts of water. For eye contact, flush immediately and for a prolonged period with copious amounts of water and consult a physician.
- In case of ingestion seek medical attention.
- Uncured monomers in this resin material may cause an allergic reaction or skin irritation in some individuals. Do not use in individuals with allergies to these materials.
- Any non-sterilizable item that is handled in the delivery of the dental service should be disinfected by standard dental office hygiene procedures.
- Do not allow saliva or water contamination of unset or set material, or of any step in the subsequent restorative procedure after applying SeptoBond.
- SeptoBond is flammable.

### STORAGE:

- The system is designed to be stored between 35°-50°F (2°-10°C).
- Replace cap immediately after dispensing.
- Allow equilibration to room temperature prior to use.
- Do not store materials in proximity of eugenol-containing products.
- Do not expose SeptoBond to elevated temperatures.
- Do not store the materials in intense light or under wet conditions.

# Image

## MICROHYBRID COMPOSITE

### RESTORATIVE MATERIAL

Light-activated dental composite for the esthetic restoration of anterior and posterior teeth

FOR PROFESSIONAL USE ONLY

INSTRUCTIONS FOR USE

### INDICATIONS:

Image Composite is indicated for direct esthetic restoration of anterior and posterior teeth; it can be used in limited occlusal function.

### DESCRIPTION:

Image Composite is a light cured resin based dental restorative material. It comes in shades matched to the Vita® Shade Guide. The material is a high viscosity, condensable composite that will resist sticking to instruments, yet will have sufficient flow when applied to the tooth to adapt to the walls of the cavity preparation.

### FORMULATION:

Approx. 75wt% (54 vol%) silanated 0.7 micron Barium glass and 0.01 micron silica; approx 25wt% BisGMA based dimethacrylate system, <1wt% catalysts, inhibitors and pigments.

### INSTRUCTIONS:

- The tooth should be cleaned by scaling and prophylaxis with flour of pumice.
- Shade selection should be performed prior to isolation and/or preparation of the teeth. Do not stare at the tabs and teeth for more than 10-15 seconds during shade selection. Use of a Vita Shade Guide is recommended.
- Isolate the teeth. Use of the rubber dam is highly recommended.
- Follow usual procedures for tooth cavity preparation.
- Base the preparation as needed. Calcium hydroxide cement may be used for covering small pulpal exposures. Any additional dentinal coverage desired can be achieved by using a glass ionomer lining cement. Eugenol containing cements must be avoided.
- Use a dentin/enamel bonding system to bond this composite to the tooth structure. Modified bonding systems can be used to bond this composite to other kinds of materials used in restorative techniques.
- Dispense the amount of composite from the syringe needed to complete the restoration(s) onto a mixing pad. Unit dose tips allow direct placement of the composite into the cavity preparation. Protect the dispensed material from direct exposure with an opaque cover. Avoid placing under the operator light during the delivery process.
- Deliver an increment of the composite directly into the preparation. Adapt the composite to the walls of the preparation using a plastic instrument. The composite should be placed into the preparation in increments no greater than 2.5 mm. It is important to use thinner layers for darker shades. Do not allow any contamination or washing between cured layers.
- With the Plasma Arc curing lights, cure each layer twice for 3 seconds. If using a conventional curing lamp, cure each layer for 30 seconds with a (halogen) blue curing light with an output around 470-480 nanometers with a minimum of 400mW/cm<sup>2</sup> intensity. An additional 20 seconds may be used from other directions. Any wedge and/or matrix can be removed after the last layer of composite is cured.
- Adjust restoration to proper margins, contours and occlusion. Remove flash with a sharp instrument. If necessary adjust occlusal contacts and finish using a fine diamond, multifluted carbide finishing bur, or finishing discs and strips. Various polishing paste are available which can enhance the polish.
- After the polishing procedure is completed, microscopic surface defects and marginal bonding that is stressed can be resealed. The composite surface can be cleaned and the surrounding enamel re-etched with phosphoric acid, washed and dried thoroughly. The bonding agent is applied in a thin layer, dried and cured.

### WARNINGS:

- Dentists and assistants should wear gloves and protective eye wear. Patients also should wear eye protection.
- If contact with eyes occurs, flush with copious amounts of water. If contact with skin occurs, wash with soap and water.
- Uncured monomers in this resin material may cause an allergic reaction or skin irritation in some individuals. Do not use in individuals with allergies to these materials.
- Any non-sterilizable item that is handled in the delivery of the dental service, should be disinfected by standard dental office hygiene procedures.
- Do not allow saliva or water contamination of unset materials or etched surfaces.

### STORAGE:

- The system is designed to be stored between 50° - 75°F (10°-24°C).
- Replace caps immediately after dispensing.
- If stored in a refrigerator, bring to room temperature prior to use.
- Do not store materials in proximity of eugenol-containing products.
- Do not store the materials in elevated temperatures, intense light, or under wet conditions.

# Image

## FLOWABLE COMPOSITE

Light-Cure Low Viscosity Composite

Resin based dental restorative material

FOR PROFESSIONAL USE ONLY

INSTRUCTIONS FOR USE

### INDICATIONS:

Image Light-cure, Low Viscosity Flowable Composite is designed for use in Class III and V restorations, and select IV aesthetic restorations. Its low viscosity and precise syringe delivery make it flowable, easy to place and finish, and adaptable to a variety of clinical situations. Image Flowable Composite is a 41% filled (by volume) composite based on a resin system of BIS-GMA and Triethylene Glycol Dimethacrylate. The filler has particle size ranging from submicron to 8 microns. Contains 2% Sodium Fluoride (by weight).

### INSTRUCTIONS:

- The tooth should be cleaned by scaling and prophylaxis with flour of pumice.
- Shade selection should be performed prior to isolation and/or preparation of the teeth. Do not stare at the tabs and teeth for more than 10-15 seconds during shade selection. Use of a Vita® Shade Guide is recommended.
- Isolate the teeth. Use of the rubber dam is highly recommended.
- Follow usual procedures for tooth preparation.

### Basing

- A calcium hydroxide or glass ionomer base may be used. Avoid the use of eugenol containing cements. Calcium hydroxide is recommended for pulp exposures.
- Use a dentin/enamel bonding system to bond this composite to the tooth structure. Modified bonding systems can be used to bond this composite to other kinds of materials used in restorative techniques.
- Attach a syringe tip to the syringe.

- Image Flowable Composite can be applied directly to the tooth surface by exerting gentle pressure on the syringe plunger. Stop pressure on plunger before the desired amount of material is dispensed, since the flow of material from the syringe tip continues for a short time after the pressure is released.
- Layer the restoration (2 to 2.5 mm) and light cure each layer for 30 seconds with a visible light curing light.

APPLICATION INSTRUCTIONS FOR UNIVERSAL OPAQUE ONLY:

- Use Universal Opaque after bonding agent application.
- Apply one thin coat (½ mm thick) to cover the discoloured or metal area only.
- Cure with a dental curing light, for 40 seconds, with a wavelength around 470-480 nanometers and with a minimum of 400mW/cm<sup>2</sup> intensity.

### Finishing

- After curing, remove flash with a sharp instrument. Image Flowable Composite provides excellent smoothness and may simplify finishing. If needed, finish using a diamond, carborundum stone, multi-fluted carbide bur or discs. Various polishing materials are available that give excellent results.

### WARNINGS:

- Dentists and assistants should wear gloves and protective eye wear. Patients also should wear eye protection.
- If contact with eyes occurs, flush with copious amounts of water. If contact with skin occurs, wash with soap and water.
- Uncured monomers in this resin material may cause an allergic reaction or skin irritation in some individuals. Do not use in individuals with allergies to these materials.
- Any non-sterilizable item that is handled in the delivery of the dental service, should be disinfected by standard dental office hygiene procedures.
- Do not allow saliva or water contamination of unset materials or etched surfaces.

### STORAGE:

- This system is designed to be stored between 50°-75°F (10°-24°C).
- If stored in a refrigerator, bring to room temperature prior to use.
- Do not store material in the proximity of eugenol-containing materials.
- Do not store the materials in elevated temperatures, intense light or under wet conditions.

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- Vita® is a registered trademark of Vita Zahnfabrik, Bad Sackingen, Germany