DEVICES DESCRIPTION

The Dentapen®

The Dentapen® is a cordless and motorized injection system powered by a lithium battery. The Dentapen® is a non-sterile device intended to assist medical practitioners in injecting dental anesthesia.

After inserting the anesthesia cartridge into the cartridge holder, the practitioner simply presses the top part of the finger grips for a smooth and even delivery of the product. The Dentapen® is only intended for dental anesthesia injections.

The finger grip

The finger grip holds a cartridge holder that receives the anesthesia cartridge, and activates Dentapen®. It has been designed to improve ergonomics, allowing the user to have greater freedom to operate during injections.

The cartridge holder

The cartridge holder is compatible with most 1.7/1.8ml anesthesia cartridges on the market and with standard dental needles. It fits in the finger grips and can still rotate while assembled. The cartridge holder and finger grips are designed for use with the Dentapen® exclusively and therefore cannot be used with any other third party device.

PRECAUTIONS AND RECOMMENDATIONS

Precautions for use

- Please read the following instructions carefully and follow the proper safety measures in order to perform a successful injection:
- Only a qualified dentist should use the Dentapen®.
- The use of the Dentapen® and the results thereof are the exclusive responsibility of the practitioner.
- The practitioner is the only one who can evaluate the indications and contraindications relative to the technique to employ, the anesthesia to be administered and the appropriate needle to use.
- In all cases, the medical practitioner shall proceed only after having fully familiarized himself or herself with both the product and its educational material, including instructions for use.
- The cartridge holder and finger grip are for a single patient and a single session use only.
- Do not use if the package is opened or damaged.
- Do not touch the keyboard while inserting the battery cap to the Dentapen®.

Adverse reactions and contraindications

It is the practitioner's full responsibility to read and follow those Instructions for use in order to avoid any adverse reactions or complaints. Possible sides effects and contraindications can result from the use of products in conjunction with the Dentapen® (e.g anesthesia cartridge, needle, etc.).

General recommendations

- Keep the Dentapen® out of reach of children, it is for professional use only.
- The equipment cannot be used where there is an inflammable mixture of air, oxygen or nitrogen.
- Store the Dentapen® away from heat and sun exposure.
- Keep the Dentapen® away from water and chemical products.
- Always check whether the device is functioning properly before use.
- If the On/Off light turns orange or red, change the battery.
- The Dentapen® turns off automatically, when not in use, after 4 minutes.
- Do not sterilize the following components: Dentapen® and battery.
- Wearing single-use gloves is highly recommended for the use of the equipment.
- No alteration of this equipment is allowed.
- Dispose of the cartridge holder, finger grips or battery according to your local environmental regulations.
- The battery used in this device may present a fire or chemical burn hazard if mistreated. Do not open battery, dispose of in fire, put in backwards, mix with used or other battery types, short circuit, disassemble, heat above 100°C (212°F), or incinerate.
- The device emits electromagnetic radiation at levels inferior to the limits recommended by the pertinent laws and regulations in force. Please refer to the EMC information section.
- Do not use the Dentapen® near electronic equipment which emits electromagnetic fields. Portable and mobile RF communications equipment can affect the Dentapen®.
- The Dentist is the exclusive responsible to check the proper aspect of the aspiration O-ring.

ACCESSORIES

Battery cap
Button Battery Lithium 6V /170 mAh
The Dentapen® must be used with a certified battery IEC60086-4 only

Finger Grip: PEEK
Cartridge holder: PSU

CLEANING

Between uses, appropriate care should be taken to ensure that the Dentapen® is clean.

- After each use, the physician shall:
- 1) Replace the needle protection, remove the needle, the cartridge holder and finger grips, the cartridge.
 - 2) Discard the cartridge and the needle.
 - 3) Press the return button (so the lead screw goes to its original configuration)
 - 4) Wipe carefully the Dentapen® with a damp cloth and a 1:100 dilution of 5.25%-6.15% sodium hypochlorite to remove visible organic residue (e.g., residue of blood) and inorganic salts.

Warning:

- Clean the Dentapen® as soon as practical after use.
- Do not immerse or soak the Dentapen® in liquids.
- Do not use an automatic washer to clean the Dentapen®.
- Inspect Dentapen® surfaces for breaks in integrity that would impair cleaning. Do not use a Dentapen® that no longer functions as intended or cannot be properly cleaned.

HANDLING CONDITIONS

It is recommended to store and transport the Dentapen® in its original packaging.

OPERATING CONDITIONS 0 - 2000 m

Temperature: 20°C TO 25°C
Humidity: 60% RH AT 30°C non-condensing

TRANSPORT CONDITIONS

Temperature: -20°C TO 40°C
Humidity: 60% RH AT 30°C non-condensing

STORAGE CONDITIONS

Temperature: -20°C TO 40°C
Humidity: 60% RH AT 30°C non-condensing

ASSEMBLY AND USE

- 1) Insert the cartridge into the cartridge holder.
- 2) Engage the cartridge holder into the finger grips.
- 3) Connect the finger grips to the main body.
- 4) Rotate the assembly counter clockwise to lock it onto the Dentapen®.
- 5) An audible "Click" indicates the proper attachment of the finger grips to the Dentapen®.
- 6) Attach the needle onto the cartridge holder.
- 7) Rotate the cartridge holder to orientate the needle as wished.
- 8) Press the ON/OFF button. The green light indicates that the device is turned ON.
- 9) Select the delivery speed, ramp mode or intraligamentary mode.
- 10) The injection can be performed by pressing the top part of the finger grips.
- 11) The injection will automatically stop when the anaesthesia cartridge is empty.
- 12) Pressing the top part of the finger grips while the Dentapen® is in motion immediately stop the injection.
- 13) Aspiration: after having activated the Dentapen®, the dentist can rapidly click again on the activation button while remaining in the same position. Dentapen® will then generate a self-aspiration feature by pulling backwards its plunger and holding it still for a few seconds. This allows the Dentist to notice if blood gets into the cartridge. If so, the Dentists should remove the needle from this site and replace the cartridge by a new one to perform another injection.
Note: always check the plunger's O-ring for cracks or deterioration prior to initiate the aspiration function as a properly maintained O-ring is necessary for effective functioning of this feature. If the O-ring is cracked or deteriorated, replace prior to use by extending the plunger while no accessory is assembled, remove the old one and slide a new one into the O-ring groove at the end of the plunger.

The Device Control Buttons and Control Lights

- Activation Button
- Device ON/OFF button, light indicator and counter-pressure indication

• **Dentapen® inactivated:**

- Green led: battery charged
- Orange led: weak battery
- Red led: battery to be replaced

• **Dentapen® activated:**

- Intraligamentary mode:**
- Blinking blue led: low pressure
- Slow blinking purple led: expected pressure
- Fast blinking purple led: high pressure

Other than intraligamentary modes:

- Slow blinking blue led: expected injection pressure

Note: in case of excessive pressure, Dentapen stops and displays an orange steady led. Reposition the needle and initiate again the injection.


- Continuous flow mode speed setup button and light indicator
- Ramp & intraligamentary activation button and light indicator
- Original position screw button and light indicator

WARRANTY

The warranty period for the Dentapen® injector is 24 months; fingers grips and cartridge holders are covered 3 months.

LEGEND

	Manufacturer		Non-sterile
	Manufacturing date		Read usage instructions
	Lot number		Keep dry
	Serial number		Do not expose to sunlight
	Catalog number		Type BF applied part
	Caution, consult documents		No Protection against ingress of solid and liquid
	Temperature limitation (-20°C/+40°C)		Do not throw away

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EMC Information

Medical electrical equipment should be used with precautions according to EMC, and must be installed according to the EMC notices disclosed in this manual.

Table 1 – Guidance and manufacturer’s declaration – Electromagnetic emissions – for all ME equipment ad ME systems		
The Dentapen® is intended for use in the electromagnetic environment specified below. The customer or the user of the Dentapen® should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	Dentapen® uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Dentapen® is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Table 2 – Guidance and manufacturer’s declaration – Electromagnetic immunity – for all ME equipment ad ME systems			
The Dentapen® is intended for use in the electromagnetic environment specified below. The customer or the user of the Dentapen® should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. Dentapen® is a self-powered (non-rechargeable battery) handy device that is not connected to the public low-voltage power supply network.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. Dentapen® is a self-powered (non-rechargeable battery) handy device that is not connected to the public low-voltage power supply network.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Model Dentapen® requires continued operation during power mains interruptions, it is recommended that the Dentapen® be powered from an uninterruptible power supply or a battery. Dentapen® is a self-powered (non-rechargeable battery) handy device that is not connected to the public low-voltage power supply network.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 4 – Guidance and manufacturer’s declaration – Electromagnetic immunity – for ME equipment and ME systems that are not life supporting			
The Dentapen® is intended for use in the electromagnetic environment specified below. The customer or the user of the Dentapen® should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the Dentapen®, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended Separation Distance d = (3,5/V1) √P2 d = 1.17 √P2 80 MHz to 800 MHz d = 2.33 √P2 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3V/m	

Table 6 – Recommended separation distances between portable and mobile RF communications equipment and the Dentapen®		
Dentapen® is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Dentapen® can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Dentapen® – as recommended below, according to the maximum output power of the communications equipment.		
Rated Maximum Output Power of Transmitter (W)	Separation Distance According to Frequency of Transmitter (m)	
	80 MHz à 800 MHz d = 1.17 √P	800 MHz à 2,5 GHz d = 2.33 √P
0.01	0.12	0.23
0.1	0.37	0.74
1	1.17	2.33
10	3.7	7.37
100	11.70	23.30