

**Root Canal Therapy should not leave you wondering whether your patients will be in pain after they have left your office.**

## **BioRoot™ RCS** **Bioceramic Root Canal Sealer**



### **EASY**

**Single cone placement**

Let the sealer be the filling material

Retreatable

### **EFFECTIVE**

**No shrinkage – resin-free<sup>1</sup>**

Dentin structure mineralization through hydroxy-apatite formation<sup>2</sup>

Allows increase of pH value >11

**Hydrophilic** – continues the setting process in the presence of moisture

Ideal in cases of apical lesions and periapical lesions

Continues to seal at the apex after placement

### **FAST**

No time spent heating gutta percha

No time spent using lateral condensation with several cones

**Fast to mix, fast to place**

Great flowability

### **INEXPENSIVE**

Priced competitively to save you money

No need to buy additional equipment

**No need for specific coated cones**



